

TruTemp[®] by TAYLOR[®]

Instant Read 1" Dial Thermometer

INSTRUCTIONS:

1. Remove meat from oven and insert thermometer to check internal temperature. Insert at least 2" of stem into food.
2. Remove thermometer and return meat to oven for further roasting until correct temperature is reached.
3. If using a microwave, read microwave operating manual carefully for guidelines on using a thermometer.

CAUTION:

Do not leave thermometer in meat while cooking in conventional ovens or rotisseries. Their high heat will damage the plastic lens of the temperature dial. Do not allow any part of the thermometer to touch the walls or doors of microwave oven or damage may occur. Do not operate microwave oven with thermometer inside. Damage may result to both oven and thermometer.

Cleaning and care: Do not put thermometer in dishwasher. Wipe clean with a damp cloth. Clean stem with alcohol to disinfect.

ONE YEAR WARRANTY

Taylor warrants this product to be free from defects in material or workmanship, and on return, will repair or replace it at our option, without charge.



Includes a versatile
pocket clip holder

